DCOCC0047
RAC L0001512190-P4-D9-1/20 0-0-0-0
70001512190000339

R

61331761206-2

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¹ HRS/EHRA Expert Consensus on Monitoring Cardiovascular Implantable Electronic Devices (CIED), April 2008.

FIELD SAFETY NOTICE

Extended Charge Time on St. Jude Medical Ellipse™ VR/DR Implantable Cardioverter Defibrillators (ICDs)

Model Numbers: CD1277, CD1377, CD2277, CD2377 (all -36, -36Q, -36C and -36QC suffixes)

19 August, 2014

Dear Customer,

As part of St. Jude Medical's focus on reliability and safety we continuously monitor the performance of our products. This letter is intended to provide you with important information regarding our Ellipse™ family of Implantable Cardioverter Defibrillators (ICDs). During capacitor maintenance or charging for high voltage therapy, a "Capacitor Charge Time Limit reached" message may occur due to an anomaly with the high voltage capacitor. This anomaly may result in delayed delivery of high voltage therapy or prevent delivery of part or all of a programmed high voltage therapy shock. There have been no serious injuries or patient deaths reported to St. Jude Medical as a result of this anomaly. No other St. Jude Medical device models are affected.

This letter is aimed to provide you with information regarding this potential anomaly, its clinical implications and patient management recommendations.

This document contains important information for the continued safe and proper use of your equipment

Please review the following information with all members of your staff who need to be aware of the contents of this communication. It is important to understand the implications of this communication.

Please retain a copy with the equipment Instructions for Use.

The U.S. Food and Drug Administration (FDA) and Competent Authorities outside the U.S. have been notified of this issue.

If you need any further information or support concerning this issue, please contact your local St. Jude Medical Representative or Technical Support at +46 8 474 4147.

We apologize for any inconvenience this may cause and trust that this information adequately addresses any concerns you may have.

Yours Sincerely,

Mark Carlson, M.D.
Vice President, Global Clinical Affairs and Chief Medical Officer

Jeff Fecho
Vice President, Global Quality

Affected products Model Numbers Serial Numbers	<p>St. Jude Medical Ellipse™ VR/DR Implantable Cardioverter Defibrillators (ICDs)</p> <p>CD1277, CD1377, CD2277, CD2377 (all -36, -36Q, -36C and -36QC suffixes)</p> <p>All serial numbers.</p>
Anomaly identified	<p>Extended charge time may present as a result of internal damage to the capacitors used in the high voltage charging circuitry of the subject devices.</p> <p>As designed, the device will deliver the available energy on the capacitors once the charge time limit of 32 seconds is reached, even if the energy is less than the programmed value.</p> <p>This anomaly may occur during capacitor maintenance or charging for high voltage therapy, and may result in delayed delivery of high voltage therapy and/or delivery of part of a programmed high voltage therapy shock.</p>
Detection of the anomaly	<p>If this condition occurs, it is detectable via a vibratory patient alert and, for patients enrolled and actively being followed, a Merlin.net notification. Additionally, upon device interrogation with the Merlin programmer, an alert message will indicate "Capacitor charge time limit reached" on a specific date.</p> <p>Approximately 97% of Ellipse extended charge time events reported to St. Jude Medical have been detected during capacitor maintenance with the remainder detected during defibrillation threshold (DFT) testing.</p> <p>There have been no reported cases of an Ellipse device failing to deliver high voltage therapy to a patient when needed. As designed, if this occurs during charging for high voltage therapy the device will deliver a shock with the available voltage.</p>
Root cause	<p>The capacitor geometry used in Ellipse ICDs is unique and therefore only Ellipse ICDs are impacted by this capacitor anomaly. Capacitors consist of individual layers of anodes (positive plates), cathodes (negative plates) and papers (insulation) which are stacked and aligned via a "front alignment hole." The completed stack is enclosed in a capacitor case which is then filled with electrolyte and sealed. Evidence of arcing has been observed between the anode and cathode at the front alignment hole in returned Ellipse devices that exhibited this anomaly.</p> <p>As a result of implementation of a standard manufacturing process operation, Ellipse ICDs with capacitors manufactured since August 2012 have demonstrated a lower likelihood of experiencing persistent damage</p>

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